

# SAFETY FEATURE EVALUATION FORM

## DENTAL GLOVES



Date: \_\_\_\_\_ Department: \_\_\_\_\_ Occupation: \_\_\_\_\_

Product: \_\_\_\_\_ Number of times used: \_\_\_\_\_

Please **circle** the most appropriate answer for each question. Not applicable (N/A) may be used if the question does not apply to this particular product.

- |  | agree.....disagree |
|--|--------------------|
| 1. The gloves dispense easily and quickly.....   | 1 2 3 4 5 N/A      |
| 2. The gloves are not discolored upon removal from the box.....  | 1 2 3 4 5 N/A      |
| 3. The glove <b>does not</b> have visible manufacturing defects (holes, etc).....                              | 1 2 3 4 5 N/A      |
| 4. The glove is available for a wide variety of hand sizes.....  | 1 2 3 4 5 N/A      |
| 5. The size is easily determined after it has been removed from the box. (sizes are marked differently).....   | 1 2 3 4 5 N/A      |
| 6. The glove is easy to put on, even if hands are damp.....  | 1 2 3 4 5 N/A      |
| 7. The glove retains appropriate sensitivity in the fingers.....   | 1 2 3 4 5 N/A      |
| 8. The glove protects the wrist securely.....  | 1 2 3 4 5 N/A      |
| 9. The glove is hypoallergenic.....  | 1 2 3 4 5 N/A      |
| 10. No excess powder remains after removing the glove.....   | 1 2 3 4 5 N/A      |
| 11. The glove is comfortable for extended use.....   | 1 2 3 4 5 N/A      |
| 12. The glove is easily removed.....   | 1 2 3 4 5 N/A      |
| 13. The glove allows the user to manipulate objects.....   | 1 2 3 4 5 N/A      |
| 14. The glove <b>does not</b> tear when it contacts sharp edges of teeth, restorations, or ortho brackets..... | 1 2 3 4 5 N/A      |
| 15. The glove <b>does not</b> become slippery on contact with saliva.....                                      | 1 2 3 4 5 N/A      |
| 16. Over gloves are easily donned removed.....   | 1 2 3 4 5 N/A      |
| 17. The glove <b>does not</b> have an unusual taste or odor to which patients object.....                      | 1 2 3 4 5 N/A      |

Of the above questions, which three are the most important to **your** safety when using this product?

Are there other questions which you feel should be asked regarding the safety/ utility of this product?